

# AllPure Mammalian Mitochondria Isolation Kit (for Cultured Cells)

Cat. No. ABTGDE401

Storage: **MSB at -20°C** for one year, others at 2-8°C for one year

## Description

AllPure Mammalian Mitochondria Isolation Kit for Cultured Cells provides a fast and efficient method to isolate mitochondria from cultured mammalian cells. This kit provides two options for the separation of mitochondria from cytosolic components: a reagent-based method or homogenization-based method. Reagent-based method uses a mild procedure to process single or multiple samples. The isolated mitochondria is suitable for a variety of downstream applications, including protein analysis, apoptosis, signal transduction and metabolic assays.

## Kit Contents

Component	50 rxns
Mitochondria Isolation Buffer I (MIB I)	50 ml
Mitochondria Isolation Buffer II (MIB II)	500 µl
Mitochondria Isolation Buffer III (MIB III)	65 ml
Mitochondria Storage Buffer (MSB) <b>-20°C</b>	4 ml
Protease Inhibitor Cocktail, EDTA-free (100X)	Not Provided

**Prior to use, Proteinase Inhibitor Cocktail or PMSF (not provided in the kit) should be added into MIB I and II and III.**

## Procedures

### Option A: Reagent-based Method

- Harvest  $2 \times 10^7$  cells and wash the cells with 1 ml of pre-chilled PBS. Centrifuge at  $1,000 \times g$  for 3 minutes. Discard the supernatant. Repeat the wash once.
- Add 800 µl of MIB I to cell pellet. Vortex for 5 seconds, and incubate on ice for 2 minutes.
- Add 10 µl of MIB II. Vortex for 5 seconds.
- Incubate on ice for 5 minutes. Briefly vortex every minute.
- Add 800 µl of MIB III. Invert tube 5-6 times to mix (do not vortex).
- Centrifuge at  $700 \times g$ , 4°C for 10 minutes.
- Gently transfer the supernatant to a new 2 ml microcentrifuge tube and centrifuge at  $12,000 \times g$ , 4°C for 15 minutes (for higher purity, suggest to centrifuge the supernatant at  $3000 \times g$  for 15 minutes at 4°C, but this may result in lower yield).
- Gently collect the supernatant (cytoplasmic protein). The isolated cytoplasmic proteins can be used for downstream applications or stored at -80°C.
- Add 500 µl of MIB III and resuspend the pellet.
- Centrifuge at  $12,000 \times g$ , 4°C for 15 minutes.
- Gently discard the supernatant, the pellet is mitochondria, which can be stored at -80°C or processed as following.
- (Option 1) For mitochondria will be used for protein analysis, the pellet can be dissolved and lysed with protein lysis buffer. Mitochondria or mitochondria lysate can be stored at -80°C for future use.
- (Option 2) For mitochondria used for functional analysis, MSB can be added at the ratio  $\sim 40 \mu\text{l}/1 \times 10^7$  cells. Analyze within one hour after resuspension.

# AllPure Mammalian Mitochondria Isolation Kit (for Cultured Cells)

## Option B: Homogenization

1. Harvest  $2 \times 10^7$  cells and wash the cells with 1 ml of pre-chilled PBS. Centrifuge at  $1,000 \times g$  for 3 minutes. Discard the supernatant. Repeat the wash once.
2. Add 800  $\mu$ l of MIB I to cell pellet. Vortex for 5 seconds, and incubate on ice for 2 minutes.
3. Transfer the suspension to a glass homogenizer and homogenize the cells by 30-50 strokes (to check the cell lysis efficiency, stain the cells with Trypan Blue and view under a microscope. If more than 50% cells are stained, homogenization can be stopped. Under homogenization may result in lower mitochondria yield. Over homogenization may damage mitochondria ).
4. Transfer the supernatant to a new 2 ml microcentrifuge tube.
5. Following steps are the same as the steps 5-12 described in "Reagent-based Method".

## Notes

- All steps should be carried out on ice or at  $4^{\circ}\text{C}$ .
- Use fresh cultured cells for mitochondria isolation if the isolated mitochondria will be used for functional assays.