

AllPure Blood Genomic DNA Kit

Cat. No. ABTGNA024

Storage: RNase A and Proteinase K solutions at -20°C for one year; others at room temperature ($15-25^{\circ}\text{C}$) for one year

Description

AllPure Blood Genomic DNA Kit provides a simple and convenient way to isolate high quality genomic DNA from 5-250 μl of fresh or frozen blood. Whole blood is incubated with binding/lysis buffer to release DNA. DNA is bound to silica-based column.

The isolated DNA is suitable for PCR, restriction enzyme digestion and Southern blot.

- Simple and fast, red cell lysis buffer is no longer needed.
- Complete removal of contaminants and inhibitors.
- DNA yield up to 40 μg .
- Column based purification, no organic extraction or ethanol precipitation.
- Suitable for EDTA, sodium citrate and heparin-anticoagulated fresh or frozen blood in a volume of 5 to 250 μl .

Starting material

For blood with nonnucleated erythrocytes, we suggest using 5-250 μl ; others, we suggest using 5-20 μl .

Kit Contents

Component	50 rxns	200 rxns
Binding Buffer 3 (BB3)	30 ml	110 ml
Clean Buffer 3 (CB3)	6 ml	24 ml
Wash Buffer 3 (WB3)	12 ml	2 \times 22 ml
Elution Buffer (EB)	25 ml	80 ml
RNase A (20 mg/ml) -20°C	500 μl	2 \times 1 ml
Proteinase K (20 mg/ml) -20°C	1 ml	4 \times 1 ml
Genomic Spin Columns with Collection Tubes	50 each	2 \times 100 each

Procedures

Before starting, add the indicated volume of 96%-100% ethanol into the concentrated CB3 and WB3.

Components	50 rxns	200 rxns
Clean Buffer 3 (CB3)	24 ml	96 ml
Wash Buffer 3 (WB3)	48 ml	2 \times 88 ml

All centrifugation steps are carried out at room temperature.

1. Add the appropriate volume of blood, 20 μl of Proteinase K and 500 μl of BB3 into a microcentrifuge tube. Mix for 15 seconds by vortexing, and then incubate at room temperature for 10 minutes.
Optional: If RNA-free genomic DNA is required, add 20 μl of RNase A before incubation.
2. Centrifuge briefly, add all the lysate to a spin column. centrifuge at $12,000\times g$ for 1 minute, discard the flow through.
3. Add 500 μl of CB3 (check to make sure ethanol has been added), centrifuge the tube at $12,000\times g$ for 30 seconds, discard flow through.
4. Add 500 μl of WB3 (check to make sure ethanol has been added), centrifuge the tube at $12,000\times g$ for 30 seconds, discard flow through.

AllPure Blood Genomic DNA Kit

5. Repeat step 4 once.
6. Place the spin column to a collection tube. Centrifuge the empty column at 12,000×g for 2 minutes to remove any residual WB3. Air-dry the spin column at room temperature for several minutes.
7. Place the spin column in a sterile 1.5 ml microcentrifuge tube. Add 50-200 µl of Elution Buffer (for higher yield, prewarm the buffer to 60°C) or distilled water (pH >7.0) to the center of column. Incubate at room temperature for 1 minute. Centrifuge at 12,000×g for 1 minute to elute the genomic DNA (to recover more DNA, add Elution Buffer or distilled water again to perform a second elution). For long-term storage, store the purified DNA at -20°C

Notes

- It is important not to overload the column, as this can lead to significantly lower yields than expected.
- Use fresh material and avoid repeated thawing and freezing.
- Use sterile tubes and pipette tips to avoid the contamination from DNase.