

## GC Enhancer for PAGE



Cat. No. ABTGAG101-01

Storage: at -20°C for two years

### Description

GC Enhancer can be used to increase sensitivity and specificity for GC/AT-rich template or complex template.

The stock concentration is 10X, and the working concentration can be varied between 0.5X to 5X.

### Application

- Complex templates
  - GC/AT-rich templates
- (50 µl reaction volume)

Volume of GC Enhancer (µl)	Final Concentration
2.5	0.5X
5	1X
10	2X
15	3X
20	4X
25	5X