# **AllBio Nhe I-Fast Cut**



Cat. No. ABTGREN301

Storage: at -20°C for two years Concentration: 10,000 units/ml Recognition Site

5'...GCTAGC...3' 3'...CGATCG...5'

#### Description

AllBio Nhe I-Fast Cut is expressed and purified from E.coli that carries the recombinant Nhe I gene. The molecular weight is 39.1 kDa, with the recognition site at G^CTAGC. The reaction is conducted for 5-15 minutes at 37°C, and heat-inactivated at 65°C for 20 minutes. This enzyme is not sensitive to dam or dcm methylation, but sensitive to mammalian CpG methylation.

# **Enzyme Properties**

Fast digestion in 5-15 minutes with high fidelity

#### Application

Genomic DNA, plasmid DNA, PCR product

# Kit Contents

Component	ABTGREN301-01	ABTGREN301-02
AllBio Nhe I-Fast Cut	500 units	2×500 units
10×FastCut Buffer	250 μl	250 µl
10×DNA Loading Buffer	1 ml	1 ml

# Unit Definition

One unit is defined as the amount of enzyme required to digest 1 µg of  $\lambda$  DNA in 1 hour at 37°C in a total reaction volume of 50 µl.

# **Quality Control**

**Ligation and re-cutting:** After 10-fold overdigestion with AllBio Nhe I-Fast Cut, more than 95% of the DNA fragments can be ligated with T4 DNA ligase at 25°C. Of these ligated fragments, more than 95% can be recut.

**16-Hour incubation:** A 50 μl reaction containing 1 μg of DNA and 10 units of enzyme incubated for 16 hours results in the same pattern of DNA bands as a reaction incubated for 1 hour with 1 unit of enzyme.

Blue/White screening (Terminal integrity): A DNA vector is digested at a unique site within  $lacZ\alpha$  gene with a 10-fold excess of enzyme, and then ligated, transformed and plated on X-gal/IPTG plate. Successful expression of the  $\beta$ -galactosidase indicates that  $lacZ\alpha$  gene remains integrity after cloning. A blue colony represents an intact gene, and a white colony represents an interrupted gene. To be Blue/White certified, enzymes must produce fewer than 3% white colonies.

**Exonuclease activity:** After incubation for 4 hours at 37°C, a 50  $\mu$ l reaction containing 100 units of enzyme and 1  $\mu$ g <sup>3</sup>H DNA releases less than 0.1% radioactive substance.

Endonuclease activity: After incubation for 4 hours at 37°C, a 50 μl reaction containing 15 units of enzyme with 1 μg pBR322 RFI DNA results in less than 10% conversion from RFI to RFII.

#### Storage Buffer

20~mM Tris-HCl pH7.4, 250~mM NaCl, 0.1~mM EDTA, 1.5~mM DTT,  $400~\mu\text{g/ml}$  BSA, 50% Glycerol

#### 10×FastCut Buffer

500 mM Tris-Ac pH7.9, 1 M KAc, 120 mM MgAc, 1 mg/ml BSA

www.allbioscience.com Tel: 0800-013-000 E-maill: service@allbioscience.com

# **AllBio Nhe I-Fast Cut**



# Reaction Components

Component	20 μl Volume	50 μl Volume
DNA	≤1 μg	≤2 μg
10×FastCut Buffer	2 μl	5 μl
AllBio Nhe I-Fast Cut	0.5 µl	1 μ1
$ddH_2O$	to 20 μl	to 50 μl

For digestion of  $\geq$ 2  $\mu g$  DNA , please adjust the amount of each component based on above table.

# **Reaction Condition**

Incubation for 5-10 minutes at 37°C. Enzyme is inactivated by adding 10×DNA Loading Buffer to a final concentration at 1×, or by heating at 80°C for 20 minutes.

#### Notes

- · Thaw the 10×FlyCut Buffer completely and mix well before use.
- · Low ionic strength, high enzyme concentration, glycerol concentration > 5%, or pH > 8.0 may result in star activity.

FOR RESEARCH USE ONLY

E-maill: service@allbioscience.com